

## DSPE-PEG400 carboxylic acid sodium

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	DSPE-PEG400 carboxylic acid sodium is a lipid formed by the conjugation of DSPE-PEG400 units with carboxyl groups. [DSPE-PEG400 carboxylic acid sodium] can be used in the preparation of lipid nanoparticles for drug delivery.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481