

Bovine brain Calcineurin

Chemical Properties

CAS No. : 90371-51-0

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Bovine brain Calcineurin is a serine/threonine protein phosphatase involved in numerous cellular processes and Ca ²⁺ -dependent signaling pathways.
Targets(IC50)	Calcium Channel

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481