

TMVP1446

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	TMVP1446, a peptide targeting VEGFR-3, can be labeled with cyanine 7 for use as a radiotracer in imaging metastatic sentinel lymph nodes (m-SLN) in cancer.
Targets(IC50)	VEGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481