

## p53-hDM2 cyclic peptide inhibitor 16e

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	p53-hDM2 cyclic peptide inhibitor 16e is a cell-permeable inhibitor of the p53-hDM2 interaction with an IC50 of 9.24 nM for hDM2. It is applicable in research related to cancers such as breast cancer, colorectal cancer, and cervical cancer.
Targets(IC50)	MDM-2/p53

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481