

## T4 endonuclease VII

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	T4endonucleaseVII exhibits broad substrate specificity, enabling it to bind and cleave both DNA duplexes and single-base mismatches in three- and four-way branched DNA structures.
Targets(IC50)	Endonuclease

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481