

cis,cis,cis-10,13,16-Docosatrienoyl-CoA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |   |
|---------------|---|
| Description   | cis,cis,cis-10,13,16-Docosatrienoyl-CoA (cis,cis,cis-10,13,16-Docosatrienoyl-coenzyme A) is a derivative of coenzyme A. |
| Targets(IC50) | Others  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481