

endo- α -1,5-Arabinanase, Bacteroides thetaiotaomicron

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	endo- α -1,5-Arabinanase, Bacteroides thetaiotaomicron (EC.3.2.1.99) is a glycoside hydrolase involved in the breakdown of arabinans. The enzyme exo- α -1,5-Arabinanase can cleave the main chain of arabinan polysaccharides.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481