

C18:0 DG-PEG2K-triGalNAc

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	C18:0 DG-PEG2K-triGalNAc is a functionalized PEG lipid composed of a C18:0 DG unit conjugated with a triGalNAc unit. This compound can be utilized in the preparation of lipid nanoparticles for drug delivery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481