

Secreted Phospholipase A2-IIA, Human

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Secreted Phospholipase A2-IIA, Human (EC 3.1.1.4), a member of the antimicrobial peptide family, functions by hydrolyzing phospholipids, thereby releasing lysophospholipids and fatty acids. It is particularly proficient at hydrolyzing bacterial membranes.
Targets(IC50)	Antibacterial

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481