

GEP12

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	GEP12 is a tri-receptor agonist peptide targeting GLP-1R, Y1-R, and Y2-R, with an EC50 of 17.3 nM as a GLP-1R agonist and an IC50 of 19.2 nM for binding to the receptor. GEP12 can aid in blood glucose management and weight reduction.
Targets(IC50)	Glucagon Receptor,Neuropeptide Y Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481