

GYGGGP(GPP)5GFOGER(GPP)5GPC

## Chemical Properties

CAS No. : 3063714-39-3

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	GYGGGP(GPP)5GFOGER(GPP)5GPC is a collagen type I mimic capable of binding to integrins and modifying hydrogels.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481