

JWJ-01-378

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	JWJ-01-378 is a selective TRK PROTAC degrader. It effectively degrades wild-type TRK and TPM3-TRKA fusion proteins while inhibiting downstream pERK signaling, but shows limited degradation effects on TRK inhibitor-resistant mutants and ALK fusion proteins. Additionally, JWJ-01-378 effectively suppresses cancer cell proliferation.
Targets(IC50)	ERK,PROTACs,Trk receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481