

CDK12-Cyclin K Ligand-Linker Conjugates 1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	CDK12-Cyclin K Ligand-Linker Conjugates 1 is a Target Protein Ligand-Linker Conjugate that includes a CDK12-Cyclin K ligand and a PROTAC linker, which can recruit E3 ligase. CDK12-Cyclin K Ligand-Linker Conjugates 1 is applicable for the synthesis of PROTACPP-C8.
Targets(IC50)	CDK,Target Protein Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481