

## FITC-LC-TAT (47-57) acetate

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from moisture**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	FITC-LC-TAT (47-57) acetate is an FITC-labeled TAT peptide known as a cell-penetrating peptide (CPP). It enhances the yield and solubility of heterologous proteins.
Targets(IC50)	HIV Protease

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481