

IdeS protease

Chemical Properties

CAS No. : 2764750-38-9

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	IdeS protease is an IgG-specific endopeptidase with high specificity for the lower hinge region of IgG. It inhibits the innate immune response, thereby facilitating the survival of streptococci in inflammatory environments.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481