

CSP-2503

## Chemical Properties


CAS No. : 581813-10-7

Formula: C<sub>22</sub>H<sub>26</sub>N<sub>4</sub>O<sub>2</sub>

Molecular Weight: 378.47

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CSP-2503 is a 5-HT <sub>1A</sub> receptor agonist.
Targets(IC <sub>50</sub> )	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6422 mL	13.2111 mL	26.4222 mL
5 mM	0.5284 mL	2.6422 mL	5.2844 mL
10 mM	0.2642 mL	1.3211 mL	2.6422 mL
50 mM	0.0528 mL	0.2642 mL	0.5284 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481