


Quinpirole

Chemical Properties

CAS No. :	80373-22-4
Formula:	C ₁₃ H ₂₁ N ₃
Molecular Weight:	219.33
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Quinpirole is a dopamine D2/D3 receptor agonist.
Targets(IC50)	Others,Dopamine Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5593 mL	22.7967 mL	45.5934 mL
5 mM	0.9119 mL	4.5593 mL	9.1187 mL
10 mM	0.4559 mL	2.2797 mL	4.5593 mL
50 mM	0.0912 mL	0.4559 mL	0.9119 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481