

Compound N064-0045

## Chemical Properties


CAS No. :

Formula: C<sub>38</sub>H<sub>52</sub>O<sub>3</sub>

Molecular Weight: 556.831

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N064-0045 is a useful organic compound for research related to life sciences and the catalog number is T131204.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7959 mL	8.9794 mL	17.9588 mL
5 mM	0.3592 mL	1.7959 mL	3.5918 mL
10 mM	0.1796 mL	0.8979 mL	1.7959 mL
50 mM	0.0359 mL	0.1796 mL	0.3592 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481