

4E1rcat analogue 6

Chemical Properties

CAS No. :


Formula: C₃₀H₂₀F₃N₄O₄

Molecular Weight: 515.48

Store at low temperature

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	The 4E1rcat analogue 6 has an (E)-3-(2-furanmethyl methylene)-2-pyrrolidone core structure. Its analogue 4E1rcat blocked the interaction between ATG5-ATG16L1 and ATG5-TECAIR in the in vitro binding assay, with IC ₅₀ values of 1.1 μM and 1.7 μM respectively.
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9399 mL	9.6997 mL	19.3994 mL
5 mM	0.388 mL	1.9399 mL	3.8799 mL
10 mM	0.194 mL	0.970 mL	1.9399 mL
50 mM	0.0388 mL	0.194 mL	0.388 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Xiang H, et al. Discovery of Small-Molecule Autophagy Inhibitors by Disrupting the Protein-Protein Interactions Involving Autophagy-Related 5. J Med Chem. 2023 Feb 23;66(4):2457-2476.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481