

4E1rcat analogue 7

Chemical Properties

CAS No. :


Formula: C₃₀H₂₃N₅O₅

Molecular Weight: 477.51

Store at low temperature

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4E1rcat analogue 7 possesses an (E)-3-(2-furylmethylene)-2-pyrrolidone core structure. Its analog 4E1rcat blocks the interactions between ATG5-ATG16L1 and ATG5-TECAIR in in vitro binding assays, with IC ₅₀ values of 1.1 μM and 1.7 μM, respectively.
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0942 mL	10.471 mL	20.942 mL
5 mM	0.4188 mL	2.0942 mL	4.1884 mL
10 mM	0.2094 mL	1.0471 mL	2.0942 mL
50 mM	0.0419 mL	0.2094 mL	0.4188 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Xiang H, et al. Discovery of Small-Molecule Autophagy Inhibitors by Disrupting the Protein-Protein Interactions Involving Autophagy-Related 5. *J Med Chem.* 2023 Feb 23;66(4):2457-2476.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481