


## Lucenin II

## Chemical Properties

CAS No. :	29428-58-8
Formula:	C <sub>27</sub> H <sub>30</sub> O <sub>16</sub>
Molecular Weight:	610.521
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Lucenin II is a natural product that can be used as a reference standard.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6379 mL	8.1897 mL	16.3795 mL
5 mM	0.3276 mL	1.6379 mL	3.2759 mL
10 mM	0.1638 mL	0.819 mL	1.6379 mL
50 mM	0.0328 mL	0.1638 mL	0.3276 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481