


WAY-661957

## Chemical Properties

CAS No. : 899902-32-0  
 Formula: C<sub>15</sub>H<sub>18</sub>ClNO<sub>3</sub>S  
 Molecular Weight: 327.83  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	WAY-661957 is a high-purity biochemical reagent that can be used in biological chemistry experiments and drug synthesis research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0504 mL	15.2518 mL	30.5036 mL
5 mM	0.6101 mL	3.0504 mL	6.1007 mL
10 mM	0.305 mL	1.5252 mL	3.0504 mL
50 mM	0.061 mL	0.305 mL	0.6101 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481