

Anti-Sarilumab Antibody (2B928)

Product Details

Ig Type:	IgG1, κ
Conjugation:	Unconjugated
Clone:	2B928
Purification:	Protein A affinity column

Applications

Application:	ELISA
Recommended	ELISA: 0.01-0.1 μ g/ml

Properties

Stability & Storage:	The lyophilized product store at -20°C or -80°C for 12 months. The reconstituted antibody can be stored at 2°C - 8°C for 2-3 weeks or at -20°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Sarilumab
Synonyms:	Kevzara

Research Background

Sarilumab with the trade name of Kevzara, is a human monoclonal IgG1 antibody. It binds to membrane bound and soluble interleukin 6 (IL-6) receptor forms. It was FDA-approved in for the treatment of Rheumatoid Arthritis (RA) in combination with methotrexate.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481