

Anti-Zolbetuximab Antibody (6B372)

Product Details

Ig Type:	IgG1, κ
Conjugation:	Unconjugated
Clone:	6B372
Purification:	Protein A affinity column

Applications

Application:	ELISA,Sandwich ELISA,Competitive ELISA
Recommended	ELISA: 0.01-1 $\mu\text{g/ml}$; Sandwich ELISA: 0.5-2 $\mu\text{g/ml}$; Competitive ELISA: 1-10 $\mu\text{g/ml}$

Properties

Stability & Storage:	The lyophilized product store at -20°C or -80°C for 12 months. The reconstituted antibody can be stored at 2°C - 8°C for 2-3 weeks or at -20°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Zolbetuximab
Synonyms:	anti-CLDN18;anti-Claudin18.2;anti-SFTPJ;anti-human Claudin-18.2;GC-182;anti-SFTA5;IMAB-362;Claudin-18.2 Inhibitor;claudiximab;IMAB362

Research Background

Zolbetuximab, also known by its development code IMAB362, is a novel monoclonal antibody targeting Claudin-18 isoform 2 (Claudin-18.2). Currently, it is being studied for its potential in treating gastrointestinal adenocarcinomas and pancreatic tumors.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481