

Anti-His tag Antibody-HRP (3T670)

Product Details

Ig Type:	IgG1, κ
Conjugation:	HRP
Clone:	3T670
Purification:	Protein A affinity column

Applications

Verified Activity:	A recombinant His-tagged protein was detected by Western blotting using Anti-His tag Antibody-HRP (3T670) (TMAS-00127H) at 1 μ g/mL concentration.
Application:	ELISA,WB
Recommended	ELISA: 0.05-0.5 μ g/ml; WB: 0.1-1.0 μ g/ml

Properties

Stability & Storage:	The lyophilized product store at -20°C or -80°C for 12 months. The reconstituted antibody can be stored at 2°C-8°C for 2-3 weeks. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: HHHHHH
Synonyms:	HisTag

Research Background

0

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481