

Anti-His tag Antibody-Biotin (3T670)

Product Details

| | |
|---------------|---------------------------|
| Ig Type: | IgG1, κ |
| Conjugation: | Biotin |
| Clone: | 3T670 |
| Purification: | Protein A affinity column |

Applications

| | |
|--------------------|---|
| Verified Activity: | Western blot analysis of His tagged fusion protein using Anti-His tag Antibody-Biotin (3T670) (TMAS-00127B). Predicted Size: 52 KD Observed Size: 52 KD |
| Application: | ELISA,WB |
| Recommended | ELISA: 0.02-0.2 $\mu\text{g/ml}$; WB: 0.5-1.0 $\mu\text{g/ml}$ |

Properties

| | |
|----------------------|--|
| Stability & Storage: | The lyophilized product store at -20°C or -80°C for 12 months. The reconstituted antibody can be stored at 2°C - 8°C for 2-3 weeks or at -20°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------|--|
| Immunogen: | KLH conjugated synthetic peptide: HHHHHH |
| Synonyms: | HisTag |

Research Background

0

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481