

Anti-His tag Antibody-iFluor 555 (3T670)

Product Details

Ig Type:	IgG1, κ
Conjugation:	iFluor 555
Clone:	3T670
Purification:	Protein A affinity column

Applications

Verified Activity:	Immunocytochemistry/Immunofluorescence analysis of non-transfected CHO cells (a) and His-tagged protein transfected CHO cells (b) using Anti-His tag Antibody-iFluor 555 (3T670) (TMAS-00127F555).
Application:	ICC/IF
Recommended	ICC/IF: 1-4 μ g/ml

Properties

Stability & Storage:	The lyophilized product store at -20°C or -80°C for 12 months. The reconstituted antibody can be stored at 2°C-8°C for 2-3 weeks or at -20°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: HHHHHH
Synonyms:	HisTag

Research Background

The His tag is a short amino acid sequence consists of Histidine (His) residues appended to recombinant proteins. Ni-charged resins are used for the purification of His-tagged proteins from eukaryotic and prokaryotic expression systems.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481