

## Anti-His tag Antibody-PE (3T670)

### Product Details

Ig Type:	IgG1, $\kappa$
Conjugation:	PE
Clone:	3T670
Purification:	Protein A affinity column

### Applications

Application:	FCM,ICC/IF
Recommended	FCM: 1:125-1:500; ICC/IF: 1:125-1:500

### Properties

Purity:	0.5 mg/mL
Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	KLH conjugated synthetic peptide: HHHHHH
Synonyms:	HisTag

### Research Background

The His tag is a short amino acid sequence consists of Histidine (His) residues appended to recombinant proteins.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481