

Anti-HA Tag Antibody-HRP (6X658)

Product Details

Ig Type:	IgG1, κ
Conjugation:	HRP
Clone:	6X658
Purification:	Protein A affinity column

Applications

Verified Activity:	Western blot analysis of Multiple Tag using Anti-HA Tag Antibody-HRP (6X658) (TMAS-00147H). Predicted Size: 52 kD Observed Size: 52 kD
Application:	ELISA,WB
Recommended	ELISA: 0.1-0.5 μ g/ml; WB: 0.1-0.5 μ g/ml

Properties

Stability & Storage:	The lyophilized product store at -20°C or -80°C for 12 months. The reconstituted antibody can be stored at 2°C-8°C for 2-3 weeks or at -20°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: YPYDVPDYA
Synonyms:	HATag

Research Background

0

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481