

## Anti-HA Tag Antibody-FITC (6X658)

### Product Details

Ig Type:	IgG1, $\kappa$
Conjugation:	FITC
Clone:	6X658
Purification:	Protein A affinity column

### Applications

Verified Activity:	Flow cytometric analysis of transfected HEK-293 cell expressing HA-tagged protein using Anti-HA Tag Antibody-FITC (6X658) (TMAS-00147F) or with a negative control antibody.
Application:	FCM, ICC/IF
Recommended	FCM: 0.5-2 $\mu$ g for $1 \times 10^6$ cells; ICC/IF: 0.5-2 $\mu$ g/ml

### Properties

Stability & Storage:	The lyophilized product store at $-20^{\circ}\text{C}$ or $-80^{\circ}\text{C}$ for 12 months. The reconstituted antibody can be stored at $2^{\circ}\text{C}$ - $8^{\circ}\text{C}$ for 2-3 weeks or at $-20^{\circ}\text{C}$ for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	KLH conjugated synthetic peptide: YPYDVPDYA
Synonyms:	HATag

### Research Background

0

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481