

## Anti-Phospho-Caveolin 1 (Tyr14) Antibody (5F839)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5F839
Purification:	Affinity-chromatography

### Applications

Application:	WB
Recommended	WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Caveolin 1 around the phosphorylation site of Y14
Synonyms:	LCCNS;Caveolin 1 (p-Y14);Caveolin1;p-Caveolin 1 (Tyr14);VIP21;CGL3;PPH3;BSCL3;CAV1; Caveolin 1 (p-Tyr14);p-Caveolin 1 (Y14);CAV

### Research Background

May act as a scaffolding protein within caveolar membranes. Interacts directly with G-protein alpha subunits and can functionally regulate their activity.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481