

Anti-Phospho-RBL2 (Thr986) Antibody (2P156)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Clone:	2P156
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of Phospho-p130(T986) expression in (1) Jurkat cell lysate; (2) Jurkat cell treated with phosphatase lysate.
Application:	IP,WB
Recommended	IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human RBL2 around the phosphorylation site of T986
Synonyms:	IB4;RBL2 (p-T986);PRB2;Rb2;p-RBL2 (Thr986);Retinoblastoma-related protein 2;130 kDa retinoblastoma associated protein;RBL2 (p-Thr986);Rbl2;RBR2;p130;Retinoblastoma like 2;p-RBL2 (T986);RIC128

Research Background

Key regulator of entry into cell division. Directly involved in heterochromatin formation by maintaining overall chromatin structure and, in particular, that of constitutive heterochromatin by stabilizing histone methylation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481