

Anti-Phospho-Tau (Ser356) Antibody (7M360)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7M360
Purification:	Affinity-chromatography

Applications

Application:	WB
Recommended	WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Tau around the phosphorylation site of S356
Synonyms:	MAPTL;p-Tau (S356);MTBT2;p-Tau (Ser356);PPND;Tau (p-S356);Tau (p-Ser356);MTBT1;PHF tau;MSTD;TAU;DDPAC;Mtapt;MAPT;pTau;RNPTAU

Research Background

Promotes microtubule assembly and stability, and might be involved in the establishment and maintenance of neuronal polarity. The C-terminus binds axonal microtubules while the N-terminus binds neural plasma membrane components, suggesting that tau functions as a linker protein between both.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481