

## Anti-ACAP2 Antibody (7L508)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	7L508
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of ACAP2 expression in 293 cell lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human ACAP2
Synonyms:	Centaurin beta2;ACAP 2;Cntb2;CENTB2

### Research Background

GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6).

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481