

## Anti-ADGRE5 Antibody (6X566)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Clone:	6X566
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of CD97 expression in U937 cell lysate.
Application:	FCM,ICC/IF,IHC,IHC-Fr,IP,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human ADGRE5
Synonyms:	Seven span transmembrane protein;CD97 antigen;Adhesion G protein-coupled receptor E5; ADGRE 5;Leukocyte antigen CD97;TM7LN1

### Research Background

Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481