

Anti-ADGRE5 Antibody (2U290)

Product Details

Ig Type:	IgG
Reactivity:	Human,Rat
Conjugation:	Unconjugated
Clone:	2U290
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of CD97 expression in human tonsil cell lysate.
Application:	FCM,IHC,WB
Recommended	FCM: 1:20-1:100; IHC: 1:100-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human ADGRE5
Synonyms:	Adhesion G protein-coupled receptor E5;TM7LN1;Seven span transmembrane protein; Leukocyte antigen CD97;CD97 antigen;ADGRE 5

Research Background

Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481