

Anti-AGP 1 Antibody (1C151)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1C151
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of ORM1 expression in human testis lysate.
Application:	IHC,WB
Recommended	IHC: 1:100-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human AGP 1
Synonyms:	AGP1;Orosomuroid 1;Alpha 1 acid glycoprotein;HEL S 153w;ORM;AGP A;alpha-1-AGP

Research Background

Functions as a transport protein in the blood stream. Binds various ligands in the interior of its beta-barrel domain. Also binds synthetic drugs and influences their distribution and availability in the body. Appears to function in modulating the activity of the immune system during the acute-phase reaction.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481