

Anti-Apolipoprotein B 100 Antibody (7X684)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7X684
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of ApoB expression in human serum lysate.
Application:	FCM,IP,WB
Recommended	FCM: 1:20-1:100; IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Apolipoprotein B 100
Synonyms:	Apolipoprotein B48;ApoB 48;Apolipoprotein B 100;ApoB 100;FLDB

Research Background

Apolipoprotein B is a major protein constituent of chylomicrons (apo B-48), LDL (apo B-100) and VLDL (apo B-100). Apo B-100 functions as a recognition signal for the cellular binding and internalization of LDL particles by the apoB/E receptor.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481