

Anti-ARHGAP5 Antibody (3X130)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	3X130
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of RhoGAP expression in Hela cell lysate.
Application:	FCM,IHC,WB
Recommended	FCM: 1:20-1:100; IHC: 1:100-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human ARHGAP5
Synonyms:	p190 B;Rho GTPase activating protein 5;p190 RhoGAP;RhoGAP5;ARHGAP 5;RhoGAP

Research Background

GTPase-activating protein for Rho family members.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481