

Anti-Bcl2L14 Antibody (3L963)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	3L963
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of BCLG expression in (1) Hela cell lysate; (2) NIH/3T3 cell lysate; (3) PC12 cell lysate.
Application:	FCM,ICC/IF,IHC,IHC-Fr,IP,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Bcl 2 L14
Synonyms:	BCL2L14;testicular tissue protein Li 26;Bcl 2 L 14;BCLG;BCL2 like 14;Apoptosis facilitator Bcl 2 like protein 14;Apoptosis regulator BCL G

Research Background

Plays a role in apoptosis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481