

Anti-BRAT1 Antibody (1L946)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1L946
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of BRAT1 expression in HeLa cell lysate.
Application:	FCM,ICC/IF,IP,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human BRAT1
Synonyms:	BRAT 1;C7orf27;BAAT1;RMFSL

Research Background

Involved in DNA damage response; activates kinases ATM, SMC1A and PRKDC by modulating their phosphorylation status following ionizing radiation (IR) stress.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481