

Anti-CDCA3 Antibody (2E875)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	2E875
Purification:	Affinity-chromatography

Applications

Verified Activity:	All lanes use CDCA3 Antibody, mAb, Rabbit at 1:1000 dilution for 1 hour at room temperature.
Application:	ICC/IF,WB
Recommended	ICC/IF: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human CDCA3
Synonyms:	TOME-1;Cell division cycle-associated protein 3;Gene-rich cluster protein C8;CDCA 3;C8;Trigger of mitotic entry protein 1;GRCC8;TOME1

Research Background

F-box-like protein which is required for entry into mitosis. Acts by participating in E3 ligase complexes that mediate the ubiquitination and degradation of WEE1 kinase at G2/M phase.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481