

Anti-CDK16 Antibody (4V302)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	4V302
Purification:	Affinity-chromatography

Applications

Application:	WB
Recommended	WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human CDK16
Synonyms:	PCTGAIRE;PCTAIRE;PCTAIRE1;PCTK1;CDK 16

Research Background

May play a role in signal transduction cascades in terminally differentiated cells.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481