

## Anti-Cholinesterase Antibody (7V244)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7V244
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of BCHE expression in A375 cell lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Cholinesterase
Synonyms:	Choline esterase II;BCHE;Butyrylcholine esterase;CHE2;Cholinesterase 1;CHE1;Cholinesterase (serum) 2;E1

### Research Background

Esterase with broad substrate specificity. Contributes to the inactivation of the neurotransmitter acetylcholine. Can degrade neurotoxic organophosphate esters.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481