

Anti-c-Myc-responsive protein RCL Antibody (6G139)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6G139
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of DNPH1 expression in HeLa cell lysate.
Application:	ICC/IF,IHC,IHC-Fr,WB
Recommended	ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human c Myc responsive protein RCL
Synonyms:	DNPH1;Rcl

Research Background

Part of a nucleotide salvage pathway that eliminates epigenetically modified 5-hydroxymethyl-dCMP (hmdCMP) in a two-step process entailing deamination to cytotoxic 5-hydroxymethyl-dUMP (hmdUMP), followed by its hydrolysis into 5-hydroxymethyluracil (hmU) and 2-deoxy-D-ribose 5-phosphate (deoxyribosephosphate).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481