

Anti-Collagen XI alpha 1 Antibody (5F52)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5F52
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of COL11A1 expression in K562 cell lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Collagen XI alpha 1
Synonyms:	COLL6;STL2;collagen XI a1;collagen XI a 1;COL11A1;STL3;Collagen alpha 1

Research Background

May play an important role in fibrillogenesis by controlling lateral growth of collagen II fibrils.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481