

Anti-c-Raf Antibody (1F342)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1F342
Purification:	Affinity-chromatography

Applications

Verified Activity:	All lanes use the c Raf Antibody, mAb, Rabbit at 1:1000 dilution for 1 hour at room temperature.
Application:	ICC/IF,WB
Recommended	ICC/IF: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human c Raf
Synonyms:	Protein kinase raf 1;cRaf;kinase Raf1;RAF1;Raf proto-oncogene serine/threonine-protein kinase

Research Background

Serine/threonine-protein kinase that acts as a regulatory link between the membrane-associated Ras GTPases and the MAPK/ERK cascade, and this critical regulatory link functions as a switch determining cell fate decisions including proliferation, differentiation, apoptosis, survival and oncogenic transformation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481