

Anti-CRISPR-Cas9 Polyclonal Antibody

Product Details

Ig Type:	IgG
Conjugation:	Unconjugated
Purification:	Immunoaffinity chromatography

Applications

Application:	WB
Recommended	WB: 1:1000-1:5000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Streptococcus pyogene (S.pyogenes) Cas9 protein.
Synonyms:	CRISPR-Cas9/Csn1;CRISPR-associated endonuclease Cas9/Csn1;SpyCas9;csn1;Cas9

Research Background

CRISPR/Cas9 system is one powerful gene editing tool to specifically modify target genes. Cas9 nuclease is an RNA-guided endonuclease that can catalyze cleavage of double stranded DNA.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481