

Anti-Cytokeratin 12 Antibody (4G146)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	4G146
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of Cytokeratin 12 expression in mouse eyeball lysate.
Application:	ICC/IF,IHC,IHC-Fr,WB
Recommended	ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Cytokeratin 12
Synonyms:	KRT12;Keratin 12;K12;Cytokeratin12;type I cytoskeletal 12

Research Background

Involved in corneal epithelium organization, integrity and corneal keratin expression.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481