

Anti-DAZ associated protein 1 Antibody (8B316)

Product Details

| | |
|---------------|-------------------------|
| Ig Type: | IgG |
| Reactivity: | Human,Mouse,Rat |
| Conjugation: | Unconjugated |
| Clone: | 8B316 |
| Purification: | Affinity-chromatography |

Applications

| | |
|--------------------|---|
| Verified Activity: | Western blot analysis of DAZAP1 expression in Jurkat cell lysate. |
| Application: | WB |
| Recommended | WB: 1:1000-1:2000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------|---|
| Immunogen: | A synthesized peptide: human DAZ associated protein 1 |
| Synonyms: | Dazap1 |

Research Background

RNA-binding protein, which may be required during spermatogenesis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481